

PRODUCT SPECIFICATIONS



Source	Xe Pulsed Lamp (U-Shape)
Spectrum of emission	400 + 1200nm @PL 400 500 + 1200nm @PL 500 550 + 1200nm @PL 550
Treated area	46 mm x 10 mm (4.6 cm ²) 46 mm x 18 mm (8.3 cm ²)
Pulse duration	5 + 50 ms
Pulse frequency	1 Hz (max)
Fluence	3 + 30 J/cm ² (4.6 cm ²) 3 + 16 J/cm ² (8.3 cm ²)
Electric supply	110 Vac - 10 A - 50/60 Hz
Interface	LCD colour touch screen
Dimensions	8" (20 cm) (H), 19" (48 cm) (W), 21" (54 cm) (D)
Weight	40 lbs (18 kg)



U-Shape PATENTED TECHNOLOGY

MAIN FEATURES

- LARGE AREAS CAN BE TREATED DURING A SINGLE SESSION
- HIGH TOLERABLE AND NON INVASIVE TREATMENT
- HIGH SPEED AND EFFECTIVENESS (8.3 cm²/s)
- APPLICATION-SPECIFIC HANDPIECES
- OPTIMIZED FOR MULTIPLE APPLICATIONS
- EASY-TO-USE AND ECONOMICALLY-PRICED



www.cynosurelaser.com



Cynosure, Inc.
 10 Elizabeth Drive
 Chelmsford, MA 01824 USA
 Tel: 978-256-4200
 Toll-Free: 800-886-2966
 Fax: 978-256-6556