



## PRODUCT SPECIFICATIONS

|                     |  |
|---------------------|--|
| Type of laser       | Semi-conductor                                       |
| Wavelength          | 808 nm   |
| Power               | 6x1 W (max) (MEASURED AS CLASS I LASER)              |
| Variable aspiration | Frequency 0.1 + 5 Hz - Duty cycle 10 + 90%           |
| Power supply        | 110 Vac - 3 A (max) - 50/60 Hz                       |
| Dimensions          | 47" (120 cm) (H) - 19" (48 cm) (W) - 24" (60 cm) (D) |
| Weight              | 88 lbs (40 kg)                                       |

## MAIN APPLICATIONS

- MINIMIZE CELLULITE
- PRE AND POST-LIPOSUCTION PROCEDURES
- FACIAL SMOOTHING OF CELLULITE
- BODY SMOOTHING OF CELLULITE
- ENHANCES MICROCIRCULATION
- THERAPEUTIC MASSAGE

**CE**  
0459

[www.cynosurelaser.com](http://www.cynosurelaser.com)  
[www.laserfacial.biz](http://www.laserfacial.biz)



**Cynosure, Inc.**  
10 Elizabeth Drive  
Chelmsford, MA 01824 USA  
Tel: 978-256-4200  
Toll-Free: 800-886-2966  
Fax: 978-256-6556